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Literatur

1. AlRowis R, AlMoharib HS, AlMubarak A, Bhaskardoss J, Preethanath RS, Anil S. Oral fluid-based biomarkers in periodontal disease - part 2. Gingival crevicular fluid. *J Int Oral Health*. 2014; 6: 126-135.
2. Gupta G. Gingival crevicular fluid as a periodontal diagnostic indicator- II: Inflammatory mediators, host-response modifiers and chair side diagnostic aids. *J Med Life*. 2013; 6: 7-13.
3. Ebersole JL, Schuster JL, Stevens J, Dawson D, Kryscio RJ, Lin Y, Thomas MV, Miller CS. Patterns of salivary analytes provide diagnostic capacity for distinguishing chronic adult periodontitis from health. *J Clin Immunol* 2013; 33: 271-279.
4. Mäntylä P, Stenman M, Kinane DF, Tikanoja S, Luoto H, Salo T, Sorsa T. Gingival crevicular fluid collagenase-2 (MMP-8) test stick for chair-side monitoring of periodontitis. *J Periodontal Res*. 2003; 38: 436-439.
5. Sorsa T, Tervahartiala T, Leppilähti J, Hernandez M, Gamonal J, Tuomainen AM, Lauhio A, Pussinen PJ, Mäntylä P. Collagenase-2 (MMP-8) as a point-of-care biomarker in periodontitis and cardiovascular diseases. Therapeutic response to non-antimicrobial properties of tetracyclines. *Pharmacol Res*. 2011; 63: 108-113.
6. Gupta A, Govila V, Saini A. Proteomics - The research frontier in periodontics. *J Oral Biol Craniofac Res*. 2015; 5: 46-52.
7. Beklen A, Tüter G, Sorsa T, Hanemaaijer R, Virtanen I, Tervahartiala T, Konttinen YT. Gingival tissue and crevicular fluid co-operation in adult periodontitis. *J Dent Res*. 2006; 85: 59–63.
8. Beikler T, Peters U, Prior K, Eisenacher M. Gene expression in periodontal tissues following treatment. *BMC Med Genomics*. 2008; 1: 30.
9. Iyer RP, Patterson NL, Fields GB, Lindsey ML. The history of matrix metalloproteinases: milestones, myths, and misperceptions. *Am J Physiol Heart Circ Physiol*. 2012; 303: H919-H930.
10. Klein T, Bischoff R. Physiology and pathophysiology of matrix metalloproteinases. *Amino Acids*. 2011; 41: 271-290.
11. Gaffney J, Solomonov I, Zehorai, E, Sagi I. Multilevel regulation of matrix metalloproteinases in tissue homeostasis indicates their molecular specificity in vivo. *Matrix Biol*. 2015; 44-46: 191-199.

12. Li JY, Wang HL. Biomarkers associated with periimplant diseases. *Implant Dent.* 2014; 23: 607-611.
13. Sorsa T, Tjäderhane L, Salo T. Matrix metalloproteinases (MMPs) in oral diseases. *Oral Dis.* 2004; 10: 311-318.
14. Sapna G, Gokul S, Bagri-Manjrekar K. Matrix metalloproteinases and periodontal diseases. *Oral Dis.* 2014; 20: 538-550.
15. Rathore PK, Bhushan A, Sharnamma B, Bali S, Dutt P. Biomarkers: A Critical Update. *IOSR J Dent Med Sci* 2013; 8: 74-78.
16. Ghodpage PS, Kolte RA, Kolte AP, Gupta M. Influence of phase I periodontal therapy on levels of matrix metalloproteinase 1 and tissue inhibitor of metalloproteinase 1. *Saudi Dent J.* 2014; 26: 171-175.
17. Popat R, Bhavsar NV, Popat PR. Gingival crevicular fluid levels of matrix metalloproteinase-1 (MMP-1) and tissue Inhibitor of metalloproteinase-1 (TIMP-1) in periodontal health and disease. *Singapore Dent J.* 2014; 35: 59-64.
18. Ramseier CA, Eick S, Brönnimann C, Buser D, Brägger U, Salvi GE. Host-derived biomarkers at teeth and implants in partially edentulous patients. A 10-year retrospective study. *Clin Oral Implants Res.* 2015; DOI 10.1111/clr.12566
19. Alfant B, Shaddox LM, Tobler J, Magnusson I, Aukhil I, Walker C.. Matrix metalloproteinase levels in children with aggressive periodontitis. *J Periodontol.* 2008; 79: 819-826.
20. Pietruska M, Pietruski J, Skurska A, Bernaczyk A, Zak J, Zelazowska B, Dolińska E, Paniczko-Drezek A, Wysocka J. Assessment of aprotinin influence on periodontal clinical status and matrix metalloproteinases 1, 2 and their tissue inhibitors saliva concentrations in patients with chronic periodontitis. *Adv Med Sci.* 2009; 54: 239-246.
21. Mouzakiti E, Pepelassi E, Fanourakis G, Markopoulou C, Tseleni-Balafouta S, Vrotsos I. The effect of smoking on the mRNA expression of MMPs and and TIMP-1 in untreated chronic periodontitis patients: a cross-sectional study. *J Periodontal Res.* 2011; 46: 576-583.
22. Maeso G, Bravo M, Bascones A. Levels of metalloproteinase-2 and -9 and tissue inhibitor of matrix metalloproteinase-1 in gingival crevicular fluid of patients with periodontitis, gingivitis, and healthy gingiva. *Quintessence Int.* 2007; 38: 247-252.
23. Bătăiosu M, Taisescu CI, Pisoschi CG, Pascu EI, Țuculină MJ, Dăguci L, Dăguci C, BaniȚă IM. Effects of therapy with two combinations of antibiotics on the imbalance of MMP-2÷TIMP-2 in chronic periodontitis. *Rom J Morphol Embryol.* 2015; 56: 77-83.

24. Tervahartiala T, Pirilä E, Ceponis A, Maisi P, Salo T, Tuter G, Kallio P, Törnwall J, Srinivas R, Konttinen YT, Sorsa T. The in vivo expression of the collagenolytic matrix metalloproteinases (MMP-2, -8, -13, and -14) and matrilysin (MMP-7) in adult and localized juvenile periodontitis. *J Dent Res.* 2000; 79: 1969-1977.
25. Toyman U, Tüter G, Kurtiş B, Kıvrak E, Bozkurt Ş, Yücel AA, Serdar M. Evaluation of gingival crevicular fluid levels of tissue plasminogen activator, plasminogen activator inhibitor 2, matrix metalloproteinase-3 and interleukin 1- β in patients with different periodontal diseases. *J Periodontal Res.* 2015; 50: 44-51.
26. Letra A, Silva RM, Rylands RJ, Silveira EM, de Souza AP, Wendell SK, Garlet GP, Vieira AR. MMP3 and TIMP1 variants contribute to chronic periodontitis and may be implicated in disease progression. *J Clin Periodontol.* 2012; 39: 707-716.
27. Ding C, Chen X, Zhang PT, Huang JP, Xu Y, Chen N, Zhong LJ. Matrix metalloproteinase-3 -1171 5A/6A polymorphism (rs35068180) is associated with risk of periodontitis. *Sci Rep.* 2015; 5: 11667.
28. Emingil G, Tervahartiala T, Mäntylä P, Määttä M, Sorsa T, Atilla G. Gingival crevicular fluid matrix metalloproteinase (MMP)-7, extracellular MMP inducer, and tissue inhibitor of MMP-1 levels in periodontal disease. *J Periodontol.* 2006; 77: 2040-2050.
29. Kivelä-Rajamäki M, Maisi P, Srinivas R, Tervahartiala T, Teronen O, Husa V, Salo T, Sorsa T. Levels and molecular forms of MMP-7 (matrilysin-1) and MMP-8 (collagenase-2) in diseased human peri-implant sulcular fluid. *J Periodontal Res.* 2003; 38: 583-590.
30. Leppilahti JM, Sorsa T, Kallio MA, Tervahartiala T, Emingil G, Han B, Mäntylä P. The utility of gingival crevicular fluid matrix metalloproteinase-8 response patterns in prediction of site level clinical treatment outcome. *J Periodontol* 2015; 86: 777-787.
31. Meissen R, Mintcheva M, Netuschil L. Matrix-Metalloproteinase-8-Spiegel in der periimplantären Sulkusflüssigkeit an Titan und Zirkonnitridoberflächen. *Int J Parodontologie & Restaurative Zahnheilk* 2014; 34: 91-95.
32. Ehlers V. MMP-8-Messung bei Patienten mit chronischer Parodontitis und Schwangerschaftsgingivitis. *Dtsch Zahnärztl Z* 2008; 63: 206-208.
33. Kinney JS, Morelli T, Oh M, Braun TM, Ramseier CA, Sugai JV, Giannobile WV. Crevicular fluid biomarkers and periodontal disease progression. *J Clin Periodontol.* 2014; 41: 113-120.
34. Izadi-Borujeni S, Mayer M, Eickholz P. Activated matrix metalloproteinase-8 in saliva as diagnostic test for periodontal disease? A case-control study. *Med Microbiol Immunol.* 2015; 204: 665-672.

35. Gursoy UK, Könönen E, Huuromonen S, Tervahartiala T, Pussinen PJ, Suominen AL, Sorsa T. Salivary type I collagen degradation endproducts and related matrix metalloproteinases in periodontitis. *J Clin Periodontol.* 2013; 40: 18-25.
36. Gupta N, Gupta ND, Gupta A, Khan S, Bansal N. Role of salivary matrix metalloproteinase-8 (MMP-8) in chronic periodontitis diagnosis. *Front Med.* 2015; 9: 72-76.
37. Kim HD, Sukhbaatar M, Shin M, Ahn YB, Yoo WS. Validation of periodontitis screening model using sociodemographic, systemic, and molecular information in a Korean population. *J Periodontol.* 2014; 85: 1676-1683.
38. Kraft-Neumärker M, Lorenz K, Koch R, Hoffmann T, Mäntylä P, Sorsa T, Netuschil L. Full-mouth profile of active MMP-8 in periodontitis patients. *J Periodontal Res.* 2012; 47: 121-128.
39. Sexton WM, Lin Y, Kryscio RJ, Dawson DR 3rd, Ebersole JL, Miller CS. Salivary biomarkers of periodontal disease in response to treatment. *J Clin Periodontol.* 2011; 38: 434-441.
40. Gursoy UK, Könönen E, Pradhan-Palikhe P, Tervahartiala T, Pussinen PJ, Suominen-Taipale L, Sorsa T. Salivary MMP-8, TIMP-1, and ICTP as markers of advanced periodontitis. *J Clin Periodontol.* 2010; 37: 487-493.
41. Javed F, Ahmed HB, Saeed A, Mehmood A, Bain C. Whole salivary interleukin-6 and matrix metalloproteinase-8 levels in patients with chronic periodontitis with and without prediabetes. *J Periodontol.* 2014; 85: e130-135.
42. Rathnayake N, Akerman S, Klinge B, Lundegren N, Jansson H, Tryselius Y, Sorsa T, Gustafsson A. Salivary biomarkers of oral health: a cross-sectional study. *J Clin Periodontol.* 2013; 40: 140-147.
43. Victor DJ, Subramanian S, Gnana PP, Kolagani SP. Assessment of matrix metalloproteinases-8 and -9 in gingival crevicular fluid of smokers and non-smokers with chronic periodontitis using ELISA. *J Int Oral Health.* 2014; 6: 67-71.
44. Skurska A, Dolinska E, Pietruska M, Pietruski JK, Dymicka V, Kemonia H, Arweiler NB, Milewski R, Sculean A. Effect of nonsurgical periodontal treatment in conjunction with either systemic administration of amoxicillin and metronidazole or additional photodynamic therapy on the concentration of matrix metalloproteinases 8 and 9 in gingival crevicular fluid in patients with aggressive periodontitis. *BMC Oral Health.* 2015; 15 (63).

45. Marcaccini AM, Meschiari CA, Zuardi LR, de Sousa TS, Taba M Jr, Teofilo JM, Jacob-Ferreira AL, Tanus-Santos JE, Novaes AB Jr, Gerlach RF. Gingival crevicular fluid levels of MMP-8, MMP-9, TIMP-2, and MPO decrease after periodontal therapy. *J Clin Periodontol.* 2010; 37: 180-190.
46. Pan Y, Li D, Cai Q, Zhang W, Ma J, Wang M, Wang L. MMP-9 -1562C>T contributes to periodontitis susceptibility. *J Clin Periodontol.* 2013; 40: 125-130.
47. Gonçalves PF, Huang H, McAninley S, Alfant B, Harrison P, Aukhil I, Walker C, Shaddox LM. Periodontal treatment reduces matrix metalloproteinase levels in localized aggressive periodontitis. *J Periodontol.* 2013; 84: 1801-1808.
48. Han DH, Shin HS, Paek D, Kim HD. Gingival crevicular fluid levels of matrix metalloproteinases cross-sectionally related to periodontitis and metabolic syndrome in community Koreans. *J Clin Periodontol.* 2012; 39: 1125-1131.
49. Nagasupriya A, Rao DB, Ravikanth M, Kumar NG, Ramachandran CR, Saraswathi TR. Immunohistochemical expression of matrix metalloproteinase 13 in chronic periodontitis. *Int J Periodontics Restor Dent.* 2014; 34: e79-84.
50. Hernández M, Gamonal J, Tervahartiala T, Associations between matrix metalloproteinase-8 and -14 and myeloperoxidase in gingival crevicular fluid from subjects with progressive chronic periodontitis: a longitudinal study. *J Periodontol.* 2010; 81: 1644-1652.
51. Leppilähti JM, Ahonen MM, Hernández M, Munjal S, Netuschil L, Uitto VJ, Sorsa T, Mäntylä P. Oral rinse MMP-8 point-of-care immuno test identifies patients with strong periodontal inflammatory burden. *Oral Dis.* 2011; 17: 115-122.
52. Liljestränd JM, Gursoy UK, Hyvärinen K, Pussinen P. Combining salivary pathogen and serum antibody levels improves their diagnostic ability in detection of periodontitis. *J Periodontol.* 2014; 85: 123-131.
53. Nomura Y, Shimada Y, Hanada N, Numabe Y, Kamoi K, Sato T, Gomi K, Arai T, Inagaki K, Fukuda M, Noguchi T, Yoshie H. Salivary biomarkers for predicting the progression of chronic periodontitis. *Arch Oral Biol.* 2012; 57: 413-420.
54. Nam SH, Jung HI, Kang SM, Inaba D, Kwon HK, Kim BI. Validity of screening methods for periodontitis using salivary hemoglobin level and self-report questionnaires in people with disabilities. *J Periodontol.* 2015; 86: 536-545.
55. Sánchez GA, Miozza VA, Delgado A, Busch L. Salivary IL-1 β and PGE₂ as biomarkers of periodontal status, before and after periodontal treatment. *J Clin Periodontol.* 2013; 40:1112-1117.

56. Shoval I, Kushner JA, Sukhu B, Wood R, Kiss T, Lawrence HP, Tenenbaum HC. The relationship between mouthrinse metalloproteinase (MMP-1, 8, 13) and albumin levels with the degree of oral mucositis in allogenic stem cell transplant patients. *Bone Marrow Transplant*. 2005; 36: 33-38.
57. Totan A, Miricescu D, Parlatescu I, Mohora M, Greabu M. Possible salivary and serum biomarkers for oral lichen planus. *Biotech Histochem*. 2015; 3: 1-7.
58. Chaussain-Miller C, Fioretti F, Goldberg M, Menashi S. The role of matrix metalloproteinases (MMPs) in human caries. *J Dent Res*. 2006; 85: 22-32.
59. Hedenbjörk-Lager A, Bjoerndal L, Gustafsson A, Sorsa T, Tjäderhane L, Åkerman S, Ericson D. Caries correlates strongly to salivary levels of matrix metalloproteinase-8. *Caries Res*. 2015; 49: 1-8.
60. Rathnayake N, Gustafsson A, Norhammar A, Kjellström B, Klinge B, Rydén L, Tervahartiala T, Sorsa T. Salivary matrix metalloproteinase-8 and -9 and myeloperoxidase in relation to coronary heart and periodontal diseases: A subgroup report from the PAROKRANK study (periodontitis and its relation to coronary artery disease). *PLoS One*. 2015; 10(7): e0126370.
61. Leppilähti JM, Kallio MA, Tervahartiala T, Sorsa T, Mäntylä P. Gingival crevicular fluid matrixmetalloproteinase-8 levels predict treatment outcome among smokers with chronic periodontitis. *J Periodontol*. 2014; 85: 250-260.
62. Pozo P, Valenzuela MA, Melej C, Zaldívar M, Puente J, Martínez B, Gamonal J. Longitudinal analysis of metalloproteinases, tissue inhibitors of metalloproteinases and clinical parameters in gingival crevicular fluid from periodontitis-affected patients. *J Periodontal Res*. 2005; 40: 199-207.
63. Reinhardt RA, Stoner JA, Golub LM, Lee HM, Nummikoski PV, Sorsa T, Payne JB. Association of gingival crevicular fluid biomarkers during periodontal maintenance with subsequent progressive periodontitis. *J Periodontol*. 2010; 81: 251-259.
64. Baeza M, Garrido M, Hernández-Ríos P, Dezerega A, García-Sesnich J, Strauss F, Aitken JP, Lesaffre E, Vanbelle S, Gamonal J, Brignardello-Petersen R, Tervahartiala T, Sorsa T, Hernández M. Diagnostic accuracy for apical and chronic periodontitis biomarkers in gingival crevicular fluid: An exploratory study. *J Clin Periodontol*. 2015; Nov 11. doi: 10.1111/jcpe.12479.
65. Honibald EN, Mathew S, Padmanaban J, Sundaram E, Ramamoorthy RD. Perioceutics: Matrix metalloproteinase inhibitors as an adjunctive therapy for inflammatory periodontal disease. *J Pharm Bioallied Sci* 2012; 4: 417-421.

66. Catona J, Ryan ME. Clinical studies on the management of periodontal diseases utilizing subantimicrobial dose doxycycline (SDD). *Pharmacol Res.* 2011; 63: 114-120.
67. Golub LM, Lee HM, Stoner JA, Sorsa T, Reinhardt RA, Wolff MS, Ryan ME, Nummikoski PV, Payne JB. Subantimicrobial-dose doxycycline modulates gingival crevicular fluid biomarkers of periodontitis in postmenopausal osteopenic women. *J Periodontol* 2008; 79: 1409–1418.
68. Soory M. A role for non-antimicrobial actions of tetracyclines in combating oxidative stress in periodontal and metabolic diseases: a literature review. *Open Dent J.* 2008; 2: 5-12.
69. Cafiero C, Matarasso S. Predictive, preventive, personalised and participatory periodontology: 'the 5Ps age' has already started. *EPMA J.* 2013; 4:16.
70. Ebersole JL, Nagarajan R, Akers D, Miller CS. Targeted salivary biomarkers for discrimination of periodontal health and disease(s). *Front Cell Infect Microbiol.* 2015; 5: 62.
71. Leppilahti JM, Hernández-Ríos PA, Gamonal JA, Tervahartiala T, Brignardello-Petersen R, Mantyla P, Sorsa T, Hernández M. Matrix metalloproteinases and myeloperoxidase in gingival crevicular fluid provide site-specific diagnostic value for chronic periodontitis. *J Clin Periodontol.* 2014; 41: 348-356.
72. Syndergaard B, Al-Sabbagh M, Kryscio RJ, Xi J, Ding X, Ebersole JL, Miller CS. Salivary biomarkers associated with gingivitis and response to therapy. *J Periodontol.* 2014; 85: e295-303.
73. Eickholz P. Parodontale Diagnostik und ihre therapeutische Relevanz. *Parodontologie* 2015; 26: 295-305.